

## 4 compartments blue

**Single Use Disposable Polypropylene Compartment Tray for use in Disposable Sets, Decontamination Services, Hospitals and Clinics.**



### Product Benefits

- Lightweight, strong and durable.
- Ideal for use in single use surgical and custom packs and hospital processing.
- Clear graduations for easy measuring.
- Low cost.
- Clean and ready for sterilization.
- Can be sterilized using Autoclave, Ethylene Oxide or Hydrogen Plasma.
- Easy release.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization.

This Range comprises of single use, disposable products and was introduced by MECO to speed up packing and sterilizing procedures, reducing the need for retrieval and decontamination.

It is ideal for use in Sterile Services, Decontamination Units, Operating Rooms and Clinics due to its low cost, product design, and performance.

The product range now offers single use, disposable trays, kidney dishes, gallipots, sponge bowls, quivers, jugs and medicine measures.

### Product Specification

- **220mm x 150mm x 24mm**
- **Pack Size 160**
- **The following standards are applicable:**

— **ISO 9001:2008**

### Quality Management System

— **ISO 13485:2003**

### Medical Devices - Quality Management System

- **Manufactured from high quality Polypropylene**
- **Non-sterile but clean and ready for packing**
- **Autoclavable up to 135 ° C**



### Guide to Using Products from MECO

- Products are provided non-sterile but clean and ready for packing and sterilization.
- Keep SteriLite products stored in their polythene bags to avoid contamination.
- Products are single use and should not be re-processed after use.
- This product has a shelf life of 60 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary and secondary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 10 °C - 35 °C